



MATTEO GOLD SRL
Via S. dei Poli, 29
36022 Cassola (VI) Italy
Operating office: Via Sarson, 73
36061 Bassano del Grappa (VI)
Italy
TEL: (+39) 424-504381
E-mail info@matteogold.com

MGASF

It's a blue liquid product specifically for cleaning stainless steel balls.

Composed of phosphates, corrosion inhibitors organic acids and non-ionic surfactants, contains acid phosphoric.

HEALTH AND SAFETY

The compound, containing phosphoric acid is classified by the European Community Council as irritant. Contact with eyes and skin should be avoided.

METHOD OF USE:

Dose about 500 grams of product per vibrator or barrel and turn the balls for 5-8 minutes.

Remove any dirty foam with water, repeat the operation in the case of old balls or new.